

eDentity Sequencer Qc

01/11/24

1.1.1 Percentage of bases \geq Q30

Metric: Calculate the fraction of bases whose Phred score are greater or equal to Q30.

Illumina	Elements
0.92	0.93

1.1.2 Read Length Retention after Q30 Trimming (Percentage).

Metric: (Average length after trimming / original read length) * 100

Illumina	Elements
94.67	94.46

1.2.1 Homopolymer Accuracy

Metric: Average Q score decay per additional base in homopolymer runs (e.g AAAA, TTTT)

Illumina	Elements
0.05	-0.01

1.2.2 Demultiplexing Efficiency (Percentage)

Metric: (Number of reads successfully demultiplexed / Total number of reads)* 100

N/B: The data from Elements was missing unsuccessfully demultiplexed reads thus we could not calculate its demultiplexing efficiency

	Illumina	Element
Total Reads	885,110,280	
Demultiplexed Reads	608,699,786	361,464,384 (?)
Demultiplexing Efficiency	68.77	NA

1.2.3 PhiX Control performance

Metric: Assess sequencing quality using spiked-in Phix control

	Illumina	Element
PhiX Error Rate	0.005588	0.022586
PhiX Alignment Rate	0.2520	0.0000
PhiX Coverage Uniformity:		
- Coefficient of Variation:	0.1554	0.4593
- Percentage within $\pm 20\%$ of mean coverage:	0.9332	0.4330
Additional Information:		
- Total Reads:	885,110,280	361,464,384
- Mapped Reads:	223,075,671	482
- Mean Coverage:	12,410,202.48	26.73
- Coverage Range:	12,54,247 - 14,897,715	1 - 70

1.2.4 Duplicate Read Rate (Percentage)

Metrics: Percentage of duplicate reads in a standard non-amplified library.

Only reads mapping to Phix were used here under the assumption that Phix is a “standard non-amplified library”

but it appears Phix reads from Illumina were amplified therefore they may not be a good measure of duplication rate.

	Illumina	Element
Mapped Reads:	223,075,671	482
Duplicate reads:	219,577,678	27
Duplicate rate	0.9843	0.0560

1.2.5 GC Bias

Metric: Deviation of coverage in GC-rich regions compared to AT-rich regions

Only reads mapping to Phix were used to compute GC Bias

Illumina	Elements
3074930	6.24